

# LUBRICANT RINSING SOLUTION FOR SURGICAL INSTRUMENTS

## Anios RHW

### A NEUTRAL RINSING SOLUTION WITH LUBRICANT EFFECT FOR SURGICAL AND MEDICAL INSTRUMENTS

A combined product which ensures quick drying without trace. Anios RHW has been specially designed for use in combination with Aniosyme DLM or Aniosyme DLM Maxi.

- Lubricant with drying activator
- Usable with all qualities of water
- Highly compatible with polysulphones



# Anios RHW

A neutral rinsing solution with lubricant effect for surgical and medical instruments



Anios RHW  
5L can

## PRODUCT CHARACTERISTICS

- Colourless clear liquid
- ph pure product at 20°C: min 2 - max. 3
- ph diluted at 0.2%: about 7
- Density at 20°C: 1.03
- Conductivity at 0.2%: 180,5 µS/cm

## PACKAGING AND ORDERING INFORMATION

PRODUCT	PACKAGING UNIT	ORDER CODE
Anios RHW	2 x 5 L can	XXXXXX
Anios RHW	2 x 10 L can	XXXXXX

## APPLICATION AREA

A neutral rinsing solution with lubricant effect for surgical and medical instruments. A combined product which ensures quick drying without trace.

## INSTRUCTIONS FOR USE

Concentrated solution. Anios RHW is used diluted from 0.03 to 0.3% in warm or cold water (0.3 to 3ml per litre of water), according to the type of machine. Anios RHW is compatible with the cycles of thermic disinfection.

## INGREDIENTS

Non ionic and anionic surfactants, isopropanol, excipients.

## STORAGE

Dangerous. Respect the precautions for use (drawn up according to the European rules in force regarding the classification and labelling of chemical products). Storage : between +5°C and +25°C.

## CERTIFICATES

On demand we are pleased to provide you with certificates.

## MEDICAL DEVICE USE SAFELY

Please read label and product information.



FOR PROFESSIONAL USE ONLY.



**ANIOS**  
Pavé du Moulin  
59260 Lille-Hellemmes  
France

Distributed by:  
**ECOLAB DEUTSCHLAND GmbH**  
Ecolab-Allee 1  
40789 Monheim am Rhein  
Germany

**ECOLAB**<sup>®</sup>

www.ecolab.com